



Vista 2000

α/β Control Software for Radiochemistry

Highlights

Full QC -

Trends, Efficiency, Background

Multiple Systems -

Any mix to 64 detectors

Flexible Reporting -

Standard plus user customized (no database knowledge needed!)

Unified platform -

Same software controls single detector, multi-detector, automatic, and manual systems

Self configuring -


Easy field expansion of hardware

Mass Attenuation -

Individual calibration curves for each count routine; multiple curve fits to choose from

Vista 2000

Control any mix of systems to 64 detectors total - MDS, WPC-9550, and MPC-9300 currently supported





Mass Attenuation

Calibration Method - Phil's Cal

Name: Count Mode: ☒ Active

Background: Update Method: QC Chart: Source Name:


Efficiency/Crosstalk: Update Method: QC Chart: Source Name: Use Attenuation Curve: ☒ Common Curve: ☒

Alpha: Beta: AtoB:

Attenuation Curve Info: Curve Equation: Coefficients: A: B: C: Show: Batch Name: Created:

AtoB: Coefficients: A: B: C: Show: Batch Name: Created:


y = Attenuation Factor x = mg mass



Mass attenuation correction for calculating accurate sample activity



Quality Control



Completely self contained QC - define, view, and print from within Vista 2000



Alpha/Beta Count Results Sample Activity Report

PIC MPC 9604 - B

Address: 1

Sample ID **EAG99-1690-02**
EAG99

Repeat **1**

Batch ID **1690**

Count Routine **Gross Alpha Beta EPA 900.0 (Soil)**

Detector Volts **1,515.0**

Sample Qty	0.100 g	sd	0.000 g
Residual Wt	100.000 mg	sd	0.000 mg
Count Date	10/25/1999 11:20	Collection Date 1	
		Collection Date 2	8/25/1999 00:00
		Half Life	0.00 days
		Decay Factor	1.000
Sample Count Time	60.00 mins	Background Count Time	1,000.00 mins

	Efficiency %	Attenuation Factor	Activity Divisor	Background cpm	Gross counts	Gross cpm	Net cpm
Alpha	19.930	0.349	1.000	0.054	11	0.183	0.129
sd	0.438						
Alpha to Beta	18.060	1.540					

Alpha/Beta Count Results Summary Activity Report

Count Routine **Gross Alpha Beta EPA 900.0 (Soil)**
Batch ID **1690**

Sample ID **EAG99-1690-02**
EAG99

Sample Qty	0.100 g	Count Date	10/25/1999 11:20
Residual Wt	100.000 mg	Sample Count Time	60.00 mins

Alpha	8.383 ±	7.238 pCi/g	838.257 % of MPC
Beta	10.980 ±	4.024 pCi/g	MDC is Greater Than Limit
			365.985 % of MPC
			MDC is Greater Than Limit

Sample ID **EAG99-1690-05**
EAG99

Sample Qty	0.100 g	Count Date	10/25/1999 12:27
Residual Wt	100.000 mg	Sample Count Time	60.00 mins

Alpha	13.784 ±	8.716 pCi/g	1,378.372 % of MPC
Beta	6.430 ±	3.464 pCi/g	MDC is Greater Than Limit
			214.317 % of MPC
			MDC is Greater Than Limit

Sample ID **EAG99-1690-09**
EAG99

Sample Qty	0.100 g	Count Date	10/25/1999 14:15
Residual Wt	100.000 mg	Sample Count Time	60.00 mins

Alpha	7.302 ±	6.906 pCi/g	730.235 % of MPC
Beta	5.957 ±	3.348 pCi/g	MDC is Greater Than Limit
			198.552 % of MPC
			MDC is Greater Than Limit

Standard reports are provided as well as the option of creating your own custom reports - *without requiring database knowledge*